EN: Differences between Visiofocus-Pro and Visiofocus®

Visiofocus-Pro implements a list of new functions, created to maximize its performances and its use in a sanitary contest:

Possible setting:

- setting "**nurs**": MQCS requested and mandatory every 30 minutes (when used by nurses as well as for any intensive use, like in airports etc...).

The  button is disabled.

- setting "**doct**": if the thermometer's temperature is changing, the AQCS (Automatic Quick Calibration System) is automatically activated and allows the rapid stabilization of the thermometer to the room temperature: the display shows a countdown indicating the time required.

To make a measurement is necessary to wait the end of the countdown without touching the device or **is possible to do the manual MQCS (suggested procedure)**.

In Doct mode the  button is enabled.

- **air conditioning On or Off. Choose "on"** in case of high air conditioning in the room;





- **Temperature level alarm**: it's possible to choose the temperature threshold beyond which the thermometer alternates the projection of Hi.2 and the temperature for 7 seconds. The possible threshold levels are $\geq 37,0^{\circ}\text{C}$ ($\geq 98.6^{\circ}\text{F}$) or $\geq 38,0^{\circ}\text{C}$ ($\geq 100.4^{\circ}\text{F}$).

- "**mem**" **button** permits to reactivate the display backlight for 3 seconds, to show the last 9 measurement and to change the settings/modes;

- **additional instructions on the battery door**;

- **special lanyard** with additional brief instructions - to be worn around the neck.

If you want to change one of the available setting proceed as follow: with thermometer in stand-by mode press the “mem” button and hold it down; after about 8 seconds the visualization on the display changes showing in rotation the following combinations: °C, °F, ORAL, RECTAL, AXILLA, CORE, DOct, NURs, 37,0°C (98,6°F), 38,0°C (100.4°F), Air- yes, not.

ITEM DISPLAYED	DESCRIPTION	MEANING	HOW TO OBTAIN IT:
	on the display appears “°C”.	Celsius degrees is the chosen measurement unit	When “°C” appears on the display release the “mem” button.
	on the display appears “°F”.	Fahrenheit degrees is the chosen measurement unit	When “°F” appears on the display release the “mem” button.
	Using the “face” button on the display appears “ORAL”.	The thermometer is set in oral reference. The temperature is always measured at in the middle of the forehead. The thermometer will provide a temperature comparable to the temperature measured with a traditional thermometer in the mouth.	When “ORAL” appears on the display release the “mem” button.
	Using the “face” button on the display appears “RECTAL”.	The thermometer is set in rectal reference. The temperature is always measured at in the middle of the forehead. The thermometer will provide a temperature comparable to the temperature measured with a traditional thermometer in the rectum.	When “RECTAL” appears on the display release the “mem” button.
	Using the “face” button on the display appears “AXILLA”.	The thermometer is set in axilla reference. The temperature is always measured at in the middle of the forehead. The thermometer will provide a temperature comparable to the temperature measured with a traditional thermometer under the armpit.	When “AXILLA” appears on the display release the “mem” button.
	Using the “face” button on the display appears “CORE”.	The thermometer is set in core reference. The temperature is always measured at in the middle of the forehead. The thermometer will provide a temperature comparable to the core temperature.	When “CORE” appears on the display release the “mem” button.
	Using the “face” button. The display alternate the temperature to “Hi.2”.	The temperature detected exceeds the threshold levels chosen for “face” button ($\geq 37,0^{\circ}\text{C}$ ($\geq 98.6^{\circ}\text{F}$) or $\geq 38,0^{\circ}\text{C}$ ($\geq 100.4^{\circ}\text{F}$).	When “37°C (98.6°F) or 38,0°C (100.4°F) appears on the display release the “mem” button.
		HOME button disabled; MQCS requested and mandatory every 30 minutes (Highly recommended when used by nurses as well as for any intensive use, like in airports etc.)	When the “nurs” message appears on the display release the “mem” button.
	When you turn on the device, the display shows "doct" before the measurement	HOME button enabled; AQCS automatic and optional MQCS (suggested for doctors' use)	When the “doct” message appears on the display release the “mem” button.
	At the end of the measurement shows "air".	Air conditioning adjustment The ”Air” function (yes or not) allows to minimize the cooling effect of intense air conditioning on the subject/patient.	When the display will show “yes” release the “mem” button.. If you want disable the Air conditioning mode, release the button when “not” appears on the display.

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