



NM^{IN}

NANOMEDICINES INNOVATION NETWORK
RÉSEAU D'INNOVATION NANOMÉDECINES

Recent Advancements in CRISPR Reagents

HQP Capacity-Building Webinar | 12 January 2021 | 1:00 – 2:00 pm PST



IDT
INTEGRATED DNA TECHNOLOGIES



RECENT ADVANCEMENTS IN CRISPR REAGENTS

ADAM CHERNICK, PHD
FIELD APPLICATION MANAGER – FUNCTIONAL GENOMICS

Jan. 12, 2021




The information in this presentation is confidential and proprietary. No part of the presentation may be used, reproduced or distributed without the express written permission of Integrated DNATechnologies, Inc.

TODAY, IDT IS THE LEADING BRAND IN DNA AND RNA SYNTHESIS



- Founded in 1987 by Dr. Joseph Walder
- Largest custom oligonucleotide manufacturer worldwide
- >1500 employees in 9 locations
- >130,000 active customers
- >95% of ordered products are manufactured and shipped in less than 24 hours




THE WORLD'S LARGEST SUPPLIER OF CUSTOM NUCLEIC ACIDS



>64,000 oligos synthesized every day


Major Product Lines

- Major R&D teams (~80 employees)
- CRISPR
- NGS
- Genotyping
- qPCR
- Synthetic Biology
- RNAi
- Bioinformatics






CUTTING AND EDITING






COMMON CRISPR ENZYMES

SpCas9



"NGG" PAM site
dsDNA cut with blunt ends
Separate crRNA + tracrRNA (annealed)
20 nt spacer


AsCas12a (Cpf1)



"TTTV" PAM site
5-base staggered cuts
Single, short crRNA
21–24 nt spacer



MANY WAYS TO DELIVER CRISPR INTO CELLS





RNP DELIVERY

Step 1
Complex



10–20 minutes

Step 2
Deliver



30–60 minutes

Alt-R™ CRISPR SYSTEM: A COMPLETE WORKFLOW



Design

CRISPR-Cas reagents

HDR reagents

Analyze

Alt-R Cas9 gRNA design tool

- Predesigned guides
- Custom designs
- Design checking

Alt-R Cas9 HDR design tool

- Friendly UI
- Empirically defined design rules
- Integration with Cas9 gRNA designs

rhAmpSeq design tool

Alt-R gRNAs

- Cas9 crRNA:tracrRNA
- Fluorescently labeled tracrRNAs
- Cas9 sgRNA
- Cas12a crRNA
- Custom ordering for any gRNA (e.g. pegRNA, Cas13)

Alt-R CRISPR proteins

- WT Cas9
- HiFi Cas9
- Cas9 nickases
- A.s. Cas12a *Ultra*
- L.b. Cas12a *Ultra*
- Fluorescently labeled Cas9

Alt-R Electroporation Enhancers

Alt-R HDR Donor Oligos

- Up to 200 nt
- Modified ssODNs

Megamer ssDNA Fragments

- Up to 2000 nt
- Sequence-verified via NGS

Linear dsDNA Fragments

- Modified to reduce blunt integration
- Up to 3000 nt
- Sequence-verified via NGS

Alt-R HDR Enhancer

- V2 coming soon!

rhAmpSeq system for CRISPR

- Multiplexed amplicon sequencing

CRISPAItRations NGS analysis tool





- Cloud-hosted UI for analysis of CRISPR on- and off-target editing



Ait-R GUIDE RNAs

Alt-R gRNA Options for *S.p.* Cas9 & *A.s.* Cas12a



Guide RNAs	Cas9 guide RNAs			Cas12a guide RNAs	
Structure	Alt-R 2-part 	Alt-R 2-part XT 	Alt-R sgRNA 	Alt-R Cas12a crRNA 	
gRNA format	Alt-R CRISPR-Cas9 crRNA & tracrRNA	Alt-R CRISPR-Cas9 crRNA XT & tracrRNA	Alt-R CRISPR-Cas9 sgRNA	Alt-R CRISPR-Cas12a crRNA	
Components	crRNA tracrRNA	crRNA XT tracrRNA	sgRNA	crRNA	
Size (nt)	36 67	36 67	100	40–44 (41 nt recommended)	
Annealing required	Yes	Yes	No	No	
Stability	++	+++	++++	+++	
Applications	<ul style="list-style-type: none"> • Cas9-expressing cells • RNP in most cell types 	<ul style="list-style-type: none"> • Co-delivery with Cas9 plasmid/Cas9 mRNA • RNP under difficult experimental conditions (e.g., high nuclease environments) 		<ul style="list-style-type: none"> • KO/KI, RNP in most cell types • Cas12a-expressing cells 	

NOW AVAILABLE! CUSTOM ORDERING TOOL FOR ANY gRNAs



[GET HELP](#) [SIGN IN](#)

0 ITEMS \$0.00 USD

[PRODUCTS & SERVICES](#) [SUPPORT & EDUCATION](#) [TOOLS](#) [COMPANY](#)

Custom CRISPR entry

Select All ACTIONS: # of Items: 1 GO BULK INPUT

1 Item Name * ✖ 📘

Scale 📘

Custom Alt-R® gRNA, 2 nmol

Sequence * (5' → 3')

5' MOD 📄 INTERNAL 📄 3' MOD 📄 BASES 📄

Bases: 0 (Min:30 Max:150)

GC: % Tm: °C 🔧 DeltaG: kcal/mole TO RNA

Step 1: Enter Sequences (1 item)

Step 2: Order Custom CRISPR Essentials (0 items)

[CONTINUE >](#)

[Show CRISPR Help](#)

- Online ordering of any custom gRNA: including Cas13, pegRNA, gRNAs for novel nuclease systems
- Lengths: 30 to 150 nt
- Scales: 2, 10, 50, & 100 nmols
- Supports: RNA bases, Alt-R modifications, 2'OMe bases, and phosphorothioate linkages



CRISPR PROTEIN ENGINEERING



IDT CRISPR PROTEIN ENGINEERING



- Cas9 vs. Cas12a (Cpf1)
- Alt-R WT Cas9 vs. **Alt-R HiFi Cas9**
 - Problem with Cas9: off-target editing
 - Bacterial mutagenesis to evolve high-fidelity Cas9
- Alt-R WT Cas12a (Cpf1) vs. **Alt-R Cas12a (Cpf1) Ultra**
 - Problem with Cpf1: low on-target editing
 - Bacterial mutagenesis to evolve enhanced activity A.s. Cas12a
 - Actively working on developing mutants of *L.b.* Cas12a with increased activity
- Fluorescent CRISPR fusion proteins: coming soon
 - **Alt-R Cas9-GFP**
 - **Alt-R Cas9-mCherry**
 - **Alt-R A.s. Cas12a Ultra-GFP**

IDT-exclusive mutants

IDT FLUORESCENT CRISPR PROTEINS RETAIN HIGH ON-TARGET ACTIVITY




- *S.p.* Cas9 (Cas9 V3, Cas9-GFP, or Cas9-mCherry) was delivered at 2.0 μ M RNP into HEK293 cells by Nucleofector™ Nucleofection (Lonza) targeting sites within the HPRT gene
- *A.s.* Cas12a RNP (Cas12a *Ultra* or Cas12a-GFP) was delivered at a suboptimal dose of 50 nM RNP into HEK293 cells by nucleofection to achieve a range of editing across the tested HPRT sites



ENRICHMENT OF EDITED CELLS BY FACS USING Cas9-GFP




- Alt-R™ Cas9 V3 or Cas9-GFP was delivered by lipofection (10 nM RNP) into HEK293 cells
- After ~18 hours, cells were sorted into three populations. High: Top 20%, Medium 80–60%, Low: Bottom 60% of cells based on GFP signal








FLUORESCENTLY-LABELED tracrRNA




- Labeled tracrRNA was duplexed with Alt-R™ crRNA XTs
- crRNA:tracrRNA duplex complexed with Alt-R™ Cas9 V3
- RNP was delivered into cells using Lipofectamine™ RNAiMAX (Thermo Fisher) at 10 nM
- HEK293 cells were incubated for ~18 h at 37°C.
- Cells were washed twice with PBS and run on a flow cytometer



Cas9 V3
+ Unlabeled
tracrRNA



Cas9 V3
+ Labeled
tracrRNA






IDT SOLUTIONS TO IMPROVE HOMOLOGY-DIRECTED REPAIR



EDITING OUTCOMES



USE THE Alt-R HDR DESIGN TOOL FOR OPTIMAL DESIGN



Alt-R HDR Design Tool

Design and order homology-directed repair (HDR) donor templates and associated Cas9 guide RNAs for genome editing human, mouse, rat, zebrafish, or *C. elegans* targets.

< BACK 1. PICK TARGET(S) 2. SPECIFY DESIGN WATCH VIDEO ▶

Visualize & specify your mutation by opening the sliders on the interactive map or input fields below. Need additional assistance? Click 'Watch video' to learn more.

NCBI Transcript Accession NM_000518.4

Exon Intron guide RNA Reference Sequence SNP/MNP Insertion Deletion

5226961 5227026

5' 3'

G T T C A C C T T G C C C C A C A G G G C A G T A A C G G C A G A C T T C T C T C A G G A G T C A G A T G C A C C A T G G T G T
C A A G T G G A A C G G G G T G T C C C G T C A T T G C C G T C T G A A G A G G A G T C C T C A G T C T A C G T G G T A C C A C A
N V K G W L A T V A S K E E P T L H V

Start edit 5226994 SNP/MNP/Deletion length 0 ZOOM TO EDIT PREVIEW TRANSLATION ◻

Mutation Select tag (optional) INSERT

Homology arms Left Right Use default homology

Number of designs 5 Add silent mutations No Yes

Guide GTAACGGCAGACTTCTCCTC

< BACK DESIGN >


For *S.p.* Cas9 designs, guidelines based on empirical data are used to recommend donors that are associated with guide RNAs that are near the desired mutation and have both high on-target and off-target scores. For *S.p.* Cas9 D10A nickase designs, guidelines based on empirical data were used to recommend paired guide RNAs that flank the desired mutation, have PAM sites oriented away from the mutation, and have cut sites 37–70 nucleotides apart.

Alt-R HDR Design tool: a novel bioinformatics tool for ssDNA HDR template design

- Human, mouse, rat, zebrafish, nematode, or custom input
- gRNA selection using IDT Alt-R gRNA design tool
- Addition of silent mutations improves rates of HDR
- Supports WT and Cas9 Nickase (D10A) strategies
- Single or multiple (batch) design
- Alt-R modified donors available in output



DONOR TEMPLATE OPTIONS




- Blunt insertion by NHEJ observed for the linear dsDNA donor, but not for ssDNA templates
- Doses titrated up until toxicity was observed.

Amaxa nucleofection
 Jurkat cells – SERPINC1 locus
 – 2 µM Alt-R Cas9 RNP, 2 µM Alt-R Cas9 Electroporation Enhancer

Donor type for long HDR	Product	Size (nt or bp)	Toxicity	HDR efficiency
ssODN	Ultramer	45-200	Low	++++
ssDNA	Megamer	201-2000	Low-Med	+++
dsDNA (linear)	gBlock	125-3000	High	++
dsDNA (plasmid)	Gene	125-3000	Very high	+



ssDNA DONOR TEMPLATES



Unmodified
PS modified
Alt-R modified

ATAATCAGGTTCTGTTCCATAAACTCTGGATTGCATTCTCgaattcACATGGAAATGCCTCTGGAGTGTATTCTCACAGAAAAGAG
A*T*AATCAGGTTCTGTTCCATAAACTCTGGATTGCATTCTCgaattcACATGGAAATGCCTCTGGAGTGTATTCTCACAGAAAAG*A*G
/Alt-R/A*T*AATCAGGTTCTGTTCCATAAACTCTGGATTGCATTCTCgaattcACATGGAAATGCCTCTGGAGTGTATTCTCACAGAAAAG*A*G/Alt-R/


SOLUTIONS TO REDUCE BLUNT INTEGRATION WITH dsDNA HDR TEMPLATES



Alt-R MODIFIED dsDNA DONORS IMPROVE HDR RATES WHILE REDUCING BLUNT INTEGRATION EVENTS FOR LARGE INSERTIONS



Assessment of Alt-R modifications on long dsDNA donors



Alt-R modified dsDNA donors mediating a 300 bp, 500 bp, or 1 kb insertion at the SerpinC1 locus (with symmetrical 100 bp homology arms).


Donors delivered into HEK293 cells at 100 nM with 2 μ M Cas9 RNP.

Cells were treated with 1 μ M Alt-R HDR Enhancer V2 for 24 hrs.

Editing outcomes assessed by NGS: Targeted amplification long-read MinION™ sequencing (Oxford Nanopore Technologies).

Analyzed via an internal data analysis pipeline.


NEW: HDR ENHANCER V2 FURTHER IMPROVES HDR RATES





rhAmpSeq FOR EDITING QUANTIFICATION


EXAMPLE CRISPR EXPERIMENT WORKFLOW






rhAmpSeq WORKFLOW

Targeted rhAmp PCR 1




CRISPAItRations FOR THE ANALYSIS OF rhAmpSeq DATA



rhAmpSeq TECHNOLOGY AND Alt-R HiFi Cas9



Verification and validation Alt-R *S.p.* HiFi Cas9 Nuclease V3




rhAmpSeq system savings:

- 8 x 40 assay designs
- ~1500 individual PCRs reduced to <96
 - Master mix
 - gDNA
- Library quantification
- Full-time equivalent hours
 - Months to days



ANALYSIS TOOL OUTPUT





THANK YOU

www.idtdna.com/CRISPR

CRISPR@idtdna.com

achernick@idtdna.com



NMIN Capacity-Building Webinar

Q&A



Adam Chernick

Field Application Manager
Integrated DNA Technologies (IDT)



Nanomedicines Innovation Network (NMIN)



<https://nanomedicines.ca>



info@nanomedicines.ca



[@NMIN_NCE](https://twitter.com/NMIN_NCE)

