



## How to install PathoRAD Extension Studio Plus (PESP-R), latest Diagnostics Centre management software from Birlamedisoft, please follow the process below.

1. Download WinRAR software from <https://www.win-rar.com/>
2. Install 32-bit winrar software.

**Skip these 3,4,5,6 options if you have installed or using PathoGold extension Studio (PES) software already in your PC, otherwise follow these instructions.**

3. Download [https://www.pathogold.com/download/setup\\_pes.rar](https://www.pathogold.com/download/setup_pes.rar)
4. Goto Download folder and extract files from above setup.
5. Right click on **pathogold\_extensionstudio** application file and click **Run as Administrator**.
6. Follow installation process, remove clicks from Installation.

Skip these 3,4,5,6 options if you have installed or using PathoGold extension Studio software already in your PC

7. Download [https://www.pathogold.com/download/setup\\_pesp-r.rar](https://www.pathogold.com/download/setup_pesp-r.rar)
8. Goto Download folder and extract files from the above (rar) file. Open folder
9. Right click on **PathoRADPLUS** application file and click run as administrator.
10. Select drive D/E/F/G any drive, change drive letter to different drive other than 'C' drive.
11. Shortcut will appear on desktop.
12. Right click on **PathoRAD** shortcut, select properties, select compatibility, choose option, **run as an administrator**, do tick mark, and apply.

### **Date format Setting**

13. Type **change date** in search. This is for to change your system date format, you will see **additional setting** option, press it, select tab **Date**, in the **short date** format type this format as **dd/MM/yyyy**

### **Firewall Setting**

14. Type **firewall** in search and select **allow an app through windows firewall**. Click change settings, then allow another app, select the file '**PathoRAD\_ExtensionStudio\_plus**' click on Add button and tick mark on **Public** and **Private** both.
15. You are all SET now.
16. Double click on PathoRAD Icon.
17. Press **Activate Demo Version**.
18. Enter Username: **birlatech**, password: **bt**



## Birlamedisoft Pvt Ltd. Process for PESP-R LIS software download and installation.

---

**This is a Demo Version Software** which is good for entering up to 20 patients to test the full functions of the software. In this software. The same software can convert into working copy once you buy / get activation key from Birlamedisoft. To activate your software, first you need to do your lab registration by visiting <https://www.pathogold.com/Home/Index#registerhere>

End of Document