



Anika. Restore Active Living.™

Tactoset®
Injectable Bone Substitute

Tactoset *Injectable Bone Substitute*

Tactoset supports cell-mediated regeneration of new bone as material is resorbed^{1,2}

Enhanced flowability & injectability


- Contains hyaluronic acid (HA)
- HA enhances delivery of the calcium phosphate

Closed mixing system

- Minimizes environmental exposure and enhances ease of use

TACTOSET MATERIAL SET

1. Pre-filled mixing syringe containing powder
2. Glass syringe containing setting solution
3. Four 1 mL Luer lock syringes
4. Female-female Luer lock connector



ANIKA SMALL BONE 15GA DELIVERY CANNULAS

1. Drillable outer cannula (with removable drillable stylet)
GAUGE 15 GA (1.83 MM)
LENGTH 60 MM
2. Push rod



TACTOSET ORDERING INFORMATION

Item Code	Description
1000040	Medium bone 11Ga delivery cannula set including: 1 drillable outer cannula with removable drillable stylet; 1 inner cannula for end delivery; 1 inner cannula for side delivery; 1 push rod
4000050	Small bone 15Ga delivery cannula set including: 1 drillable outer cannula with removable drillable stylet; 1 push rod
6000041	Tactoset material set including: 1 prefilled mixing syringe containing Tactoset powder; 1 glass syringe containing Tactoset setting solution; 4 Luer lock syringes (1 mL); 1 female-female Luer lock connector

Please note that both the Tactoset material kit item code and one of the delivery cannula item codes are needed for each procedure.

For complete product information, please visit: www.anikaifu.com

1. Data on file, Anika Therapeutics, Inc.
2. Tactoset Instructions for Use. Anika Therapeutics, Inc.; AML 500-335.

Anika Therapeutics, Inc.
32 Wiggins Ave, Bedford, MA
1-888-721-1600 · www.anika.com

Anika and Tactoset are registered trademarks of Anika Therapeutics, Inc.
©2021 Anika Therapeutics, Inc. | 04/21 | AML 900-150/B

